

Choosing and Interpreting Diagnostic Testing in an Environment of Increasing Technology



SWINE
SCHWEIN
PORCS
SUINOS
猪
豬

Wendy Witbeck
LPD Technical Service, North America
IDEXX Livestock, Poultry and Dairy

Kristen Roza-Sutherland, DVM
LPD Sales Manager, North America
IDEXX Livestock, Poultry and Dairy

Felipe Navarro, DVM, MBA
Worldwide Swine Marketing Manager
IDEXX Livestock, Poultry and Dairy

Test With Confidence™



Agenda


- Review of fundamental testing methods
- Current and emerging technologies
- Strategic use of testing for increased information
- Clinical examples of concurrent test usage
- Interpretation of parallel testing
- Summary

Options for Diagnostic Testing Increasing

- Veterinarians and Diagnosticians continue to develop and sometimes struggle with the best use of available technologies to accurately diagnose animal disease
- As faster and more accurate technologies' become available it is important to understand how these fit into the diagnostic picture
- Technology and information sharing has expanded greatly with the invention and adoption of internet resources and online scientific forums
- New developments in diagnostics can be used to support and enhance fundamental diagnostic testing results when the user understands interpretation in context

Review of Fundamental Testing Methods

- Antibody Detection
 - ELISA
 - Agar Gel Immunodiffusion (AGID)
 - Hemagglutination Inhibition (HI)
- Antigen Detection
 - Conventional PCR
 - Virus Isolation
 - ELISA
- Confirmatory Testing
 - Culture
 - Virus Neutralization
 - Necropsy




Considerations

- Advantages
 - Often less expensive to run
 - Less equipment investment
- Disadvantages
 - Specificity
 - Sensitivity
 - Time to results

Current and Emerging Technologies

- Real time PCR
- High Throughput Sequencing
- Antibody testing with Infrared Fluorophore Testing
- Next?



Considerations

- Advantages
 - Increased sensitivity
 - Efficiency of testing and time to results
 - Increased availability of Quality Controls
 - More clarity in results, less subjectivity
- Disadvantages
 - Costs of new technologies may be higher
 - New instrumentation is sometimes needed
 - Increased technical difficulty to run/ more technician training needed

Enzyme Linked Immunosorbent Assay (ELISA)

- ELISA detects antibodies after infection.
- Antibodies are found in serum, oral fluids, and GI tract (mucosal).
- The indirect ELISA provides a semi quantitative range of antibody levels when used routinely over time.

Benefits

- Low cost
- Commercial assays—high repeatability and validated test performance

Drawback

- Time to detection of seroconversion in acute infection diagnosis



PCR—Nucleic Acid Detection


- Detects genetic material of pathogen
 - DNA or RNA
- Target sequence is amplified to enhance detection
- Reported as the number of cycles necessary for detection
 - Cycle threshold value (Ct-value)
 - Lower Ct = more initial target present

Benefit

- Highly sensitive and specific

Drawbacks

- Extremely narrow target(s), may miss variants that are slightly different
- Requires continuous field sample monitoring to ensure adequate primers (i.e., RNA viruses)



Test With Confidence™

IDEXX

Interpretation of Simultaneous Antibody and Antigen Testing Methods (ex. ELISA and PCR)

Antibody/ ELISA	Antigen/ PCR	Interpretation
+	+	Active viremia <i>(circulation and infection)</i>
+	-	Previous infection <i>(exposed but no evidence of viremia at sampling)</i>
-	+	Early infection
-	-	Negative <i>(no circulation or infection)</i>

Clinical Examples of Concurrent Test Usage

Overview

- Swine Farm breaks with clinical symptoms consistent with PRRS
- Veterinarian is called to the farm to investigate
- Samples taken are Oral Fluids, Blood, Feces
- Veterinarian would like guidance on how to test these to best understand the clinical picture
- History of the farm is expected to be PRRS negative

Diagnostic Plan:

- PRRS Oral Fluid ELISA on pooled pen samples
- Concurrent PRRS PCR on blood samples
- Consideration given to the farm size and expected disease prevalence to determine optimal number of samples

Clinical Examples of Concurrent Test Usage Results

Timeline	PCR Results	ELISA Results	Interpretation
Day 1	+	-	Clinically sick pigs, expected PRRS negative premises
Week 1	+	-	Acute infection, no detectable antibody response
Week 2	+	+	Ongoing infection, pigs clearing with veterinary supportive treatment
Week 5	-	+	Acute infection cleared from herd

Ongoing Usage of Diagnostics for Status Monitoring

Interpretation of Parallel Testing

	Classification	Shedding Status (PCR)	Exposure Status (ELISA)
Category 1	Positive Unstable	+	+
Category 2	Positive Stable	?	+
Category 3	Provisional Negative	-	+
Category 4	Negative	-	-

Differences Between Antibody and Nucleic Acid Detection

Antibody

- Detects exposure and immune response to target agent
- Detects maternal antibodies, vaccination, pathogen exposure, or reinfection
- High repeatability
- Technically simple
- Cost-effective for continuous health monitoring over time

Antigen/Nucleic acid

- Real time detection of disease agent in the animal/herd
- High sensitivity to identical target agent
 - May have cross-reactions
 - May not detect new variants
- Technically demanding
- Costly


Antibody and PCR assays detect different targets and should not be used as confirmatory tests for each other.

Early in infection, the target may be present with limited antibody production, but as disease progresses, more antibodies will be present and with a low prevalence of the target agent.

Use of Results

Antibodies

- Excellent screening tool
- Gives comprehensive picture of disease status over time
- Early awareness of changes in disease circulation
- Useful for herd health management
 - Disease elimination
 - Vaccination timing
 - Identifying concurrent infections



Antigen/Nucleic acid

- Investigation for a specific agent
- Point-in-time information
 - Population sensitivity may be low if low prevalence and small sample size
 - Difficult to detect target agent late in infections
- Semi-quantitative assessment of pathogen load
- Sample size and pooling influences ability to detect

Summary

Traditional diagnostics such as ELISA and newer technologies such as PCR are complimentary to one another and are not meant to replace the other

- Surveillance/ monitoring- ELISA
 - Examples: disease eradication programs and FAD surveillance
- Disease detection/circulating pathogen – PCR
 - Example: Early disease outbreak

Know your goals!

Thank you for your time and attention!



SWINE
SCHWEIN
PORCS
SUINOS
猪
豬

References

1. Antibody. Tutorvista website. <http://images.tutorvista.com/content/immune-system/antigen-antibodies-formation.jpeg>. Accessed September 3, 2015.
2. Formation of antigen-antibody complex. Wikimedia website. <https://upload.wikimedia.org/wikipedia/commons/thumb/2/2d/Antibody.svg/255px-Antibody.svg.png>. Accessed September 3, 2015.
3. Chemical response to disease. Quizlet website. http://o.quizlet.com/dL1wRJ6C9XnhDLTBptwaFA_m.png. Accessed September 3, 2015.
4. Immune response chart and graph. Murphy KM, Travers P, Walport M. Figure 1-24. *Janeway's Immunobiology*. New York: Garland Science; 2008.
5. HI assay image. Centers for Disease Control website. <http://wwwnc.cdc.gov/eid/images/09-1733-F1.jpg>. Accessed September 3, 2015.
6. HI assay chemistry image. National Institutes of Health website. http://openi.nlm.nih.gov/imgs/512/376/2854694/2854694_pone.0010176.g001.png. Accessed September 3, 2015.
7. Micrognome. How serology works. Micrognome website. <http://micrognome.priobe.net/2013/08/how-serology-works>. Published August 25, 2013. Accessed September 3, 2015.
8. Donofrio G, Sartori C, Ravanetti L, et al. Establishment of a bovine herpesvirus 4 based vector expressing a secreted form of the bovine viral diarrhoea virus structural glycoprotein E2 for immunization purposes. *BMC Biotechnology*. 2007;7:68.
9. Cameron SO, Carman WF. The use of the OraSure collection device for hepatitis virus testing in health care settings. *J Clin Virol*. 2005;34:22–28.
10. Chandra G. ELISA. IASZOOLOGY website. <http://www.iaszoology.com/wp-content/uploads/image/ELISA.jpg>. Accessed September 3, 2015.

Proprietary Rights Information

Information in these materials is subject to change without notice. Companies, names, and data used in examples are fictitious unless otherwise noted. No part of this document may be reproduced or transmitted in any form or by any means, electronic, mechanical, or otherwise, for any purpose, without the express written permission of IDEXX Laboratories, Inc. IDEXX Laboratories, Inc. may have patents or pending patent applications, trademark, copyrights, or other intellectual or industrial property rights covering this document or subject matter in this document. The furnishing of these materials does not give a license to these property rights except as expressly provided in any written license agreement from IDEXX Laboratories, Inc.

Any recommendations contained in these materials are intended to provide general guidance only. As with any diagnosis or treatment, you should use clinical discretion with each patient based on a complete evaluation of the patient, including physical presentation and complete laboratory data. With respect to any drug therapy or monitoring program, you should refer to product inserts for a complete description of dosages, indication, interactions, and cautions. For more specific information on IDEXX products and services, please refer to the appropriate operator guides, terms of sale, product inserts, and other materials.

© 2017 IDEXX Laboratories, Inc. All rights reserved.

* All ®/TM marks are owned by IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

The IDEXX Privacy Policy is available at Idexx.com.

Test With Confidence™

